### **Purpose**

This document describes IPG and EPG calibration.

### Scope

The scope of this document includes the Model 2408 and 2412 IPGs as well as the Model 4300 trial stimulator.

#### References

- EESP 0071 24-Channel IPG Functional Specification
- EESP 0091 25 Channel EPG Functional Specification
- SWSP 0112 IPG/EPG Software Requirement Specification

## **Output Channel Calibration**

The IPG and EPG output channels are calibrated during manufacturing. This calibration hat the following characteristics:

- a. There are four lookup tables with each one containing 52 elements. Each channel has two elements in each table: one for calibrating the current source and one for calibrating the current sink.
- b. Each lookup table is for a different amplitude level. The four tables correspond to the following amplitudes: 15mA, 9mA, 6mA, and 3mA.
- c. Each element in the table contains the value for setting the corresponding output cur to the corresponding calibration current within ±7.5µA.
- d. The values in this table are used to linearly interpolate the necessary amplitude setti required to produce the desired output current.

# **High-Voltage Supply Calibration**

The IPG's high-voltage supply shall be calibrated during manufacturing. This calibration s have the following characteristics:

- a. There is a lookup table containing 64 elements. The first element in the table shall correspond to the lowest supply setting while the last element shall correspond to the highest supply setting.
- b. Each element is loaded with a value that represents the supply's output voltage for the particular setting.
- c. The IPG/EPG uses this lookup table to determine the appropriate power supply setting achieve the desired output voltage.